

1

Fill the tank to 1 or 2 users' line (for both water and ice)

2

Connect air tubing from Aquilo control unit to each air attachment on the Boots

3


Connect water tubing from Aquilo control unit to each water attachment on the Boots

4

Put on the Boots and zip until highest point

aquilo

Quick User Guide



Upper screen (water pump)



1 Press Power button on right corner (metal button) in order to switch on the device

2 Turn on water pump by pressing on the far-left button for 3 seconds, then press again to start

Set timer by either pressing + or – to increase or decrease time (standard time is 25mins)

Pulse button will run the water pump for 2 minutes on, 30 seconds off, repeat

Bottom screen (air pump)



Press Play for Standard Pressure settings

Standard setting: 75mmHg, 30min and sequential compression – bottom up (Mode 1)

Changing compression Modes

- 1 Press Set button on the top left and then the + or – button on the top right. A Red Arrow will appear under the Mode (ex. M1, M2, ...)

Mode 2

Uniform compression – all air channels are filled at the same time

Mode 3

Sequential compression – top to bottom

Mode 4

sequential compression: bottom up with 1 channel at the time; releasing previous channel

Mode 5

sequential compression: bottom to top with 2 channels at a time

It is also possible to select specific air channels by clicking on them (1CH, 2CH..)

Bottom screen (air pump)



Changing compression levels

- 1 Press OK button below Set button and then + or – to increase or decrease levels. A Red Arrow will appear under the mmHg

Minimum compression level is 20mm of mercury
Maximum compression level is 150mm of mercury

Changing time

- 1 Press OK button below Set button and then + or – to increase or decrease time

Minimum time: 20mins
Maximum time: 200mins

Press and hold Lock button, if desired, to lock settings and prevent changes to settings

- 2 Press Start/Pause button at bottom left in order to launch the compression pump

