



## FORM FIT® OA WRAPAROUND

Lightweight, comfortable  
and discreet

Form Fit OA Wraparound is a plug-and-play knee brace for low active users with mild-to-moderate symptoms of unicompartmental osteoarthritis. Its simple closure design is particularly appropriate for patients with poor hand strength or dexterity.

### INDICATIONS FOR USE:


- Mild to moderate unicompartmental knee osteoarthritis
- Patients who prefer a wraparound design
- Patients with poor hand dexterity or cognitive issues
- Low-to-medium contact/impact



## FORM FIT® OA WRAPAROUND

Form Fit OA Wraparound is a cost-effective, low-profile option for patients with mild-to-moderate symptoms of unicompartmental knee osteoarthritis.

During leg extension, the Dynamic Force Strap tightens to provide pain relief in a single pull motion. The slim polycentric hinge design and breathable material produce a discreet, comfortable fit. Dot-coded color wrap closures allow for a simple donning procedure for the user.



### FEATURES AND BENEFITS

- 1 | Easy wraparound design  
Plug-and-play fitting
- 2 | Dot coded closure wraps  
Simple donning procedure
- 3 | Non-slippery rubber tab  
Easy grip while donning
- 4 | Dynamic Force Strap  
Pain relief with a single pull motion
- 5 | Low-profile design & slim hinge  
Discreet look
- 6 | Breathable lycra  
Comfortable & consistently snug fit

### ORDERING INFORMATION

Size	Thigh Circ. <sup>1</sup>	Calf Circ. <sup>2</sup>	Right Medial	Left Medial	Right Lateral	Left Lateral
XSmall	14-16" (36-41cm)	11-12.5" (28-32cm)	B-254901001	B-254901008	B-254901015	B-254901022
Small	16-18" (41-46cm)	12.5-14" (32-36cm)	B-254901002	B-254901009	B-254901016	B-254901023
Medium	18-21" (46-53cm)	14-16" (36-41cm)	B-254901003	B-254901010	B-254901017	B-254901024
Large	21-23.5" (53-60cm)	16-18" (41-46cm)	B-254901004	B-254901011	B-254901018	B-254901025
XLarge	23.5-26.5" (60-67cm)	18-20" (46-51cm)	B-254901005	B-254901012	B-254901019	B-254901026
2XLarge	26.5-29.5" (67-75cm)	20-22" (51-56cm)	B-254901006	B-254901013	B-254901020	B-254901027
3XLarge	29.5-32" (75-81cm)	22-24" (56-61cm)	B-254901007	B-254901014	B-254901021	B-254901028

1 Thigh Circumference: 6" above mid-patella

2 Calf Circumference: 6" below mid-patella

FOLLOW ÖSSUR ON



TEL (800) 233-6263  
CANADA (800) 663-5982  
WWW.OSSUR.COM

