



COMFORT LUX PRESSURE CARE MATTRESS RANGE

RG66M150: 150mm Mattress
RG66M200: 200mm Mattress
RG66M250: 250mm Mattress



RG66M250
250mm Mattress shown

FEATURES

- 3 mattress heights available
- Multi layer memory foam construction
- Gel infused memory foam top layer
- Shaped foam body for excellent mattress flex
- Quality woven mattress coverings
- Removable washable cover with zipper access
- Secondary mattress cover to minimise foam tear
- Non-slip fabric on underside of mattress cover

SPECIFICATIONS

Overall: **2030mm (L) x 1050mm (W)**
Maximum user weight: **200kg**
Mattress height / weight:
RG66M150: **150mm / 11kg**
RG66M200: **200mm / 15kg**
RG66M250: **250mm / 19kg**

Please contact **REDGUM Brand** for further information on this product.

P: 08 9248 4180

E: orders@for-de.com.au


1 Business Way, Malaga, Western Australia 6090

Part of the **For-de Group Pty Ltd**

ICE COOL COVER

Specifically designed for "hot" sleepers and those who prefer a cooler sleep experience.

The mattress cover is made from UHMWPE fiber that has an advanced moisture wicking structure which allows for excellent airflow and fast fabric drying.



MULTI LAYER Memory Foam Construction

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.
The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

REDGUM brand strives to provide the highest quality products possible that are subject to change without notice.

ASK FOR THE BRAND

redgumbrand.com.au