


NEW PRODUCT!



Tilt top with various
ANGLE SETTINGS



Over Bed Chair Table
can be used with:

- **BEDS**
- **KATIE CHAIRS**
- **LIFT CHAIRS**

OVER BED CHAIR TABLE

RG626

FEATURES

- Large working surface
- Height adjustable platform with retaining lip
- Tilt top with various angle settings
- Fixed side platform with retainer
- Width adjustable base
- Low clearance base to fit under beds
- Swivel castors
- Can be flat packed for storage or transport

SPECIFICATIONS

(All dimensions shown are nominal)

- Large Platform: **750mm (L) x 400mm (W)**
- Side Platform Depth: **180mm (L) x 400mm (W)**
- Platform Height: **700 - 980mm**
- Base Frame: **550 - 780mm (W-internal) x 115mm (H)**
- Castors: **Ø50mm plastic twin wheel locking castors**
- Product Weight: **12kg**
- Maximum Weight Loading: **10kg - evenly distributed**

Please contact **REDGUM Brand** if you require further information on this product.

P: **08 9248 4180**
E: **orders@for-de.com.au**

ASK FOR THE BRAND

redgumbrandmobility.com.au