

CPR  
FUNCTION

**4 SECTION**  
mattress platform



# COMFORT LUX MEDICAL HIGH LOW KING SINGLE BED

RG65BHLH

## FEATURES

- 4 section mattress platform
- Mattress retainer
- Adjustable backrest
- Adjustable knee break
- Hard wired handset with lock out function
- CPR function
- Quality electrics
- Battery backup
- 4 total locking castors

## SPECIFICATIONS


- Overall: **2050mm (L) x 1120mm (W) x 370mm (H)**  
(H - at lowest setting)
- Height Adjustment: **370 - 690mm** (top of mattress to floor)
- Backrest Platform: **710mm (L) x 1050mm (W)**  
**0° - 70°** (Angle range)
- Mid Platform: **380mm (L) x 1050mm (W)**
- Knee Break Platform: **300mm (L) x 1050mm (W)**  
**0° - 40°** (Angle range)
- Calf Support Platform: **500mm (L) x 1050mm (W)**
- Mid Platform: **380mm (L) x 1050mm (W)**
- Mobility Wheel Diameter: **Ø75mm total lock brakes**
- Product Weight: **116kg**
- Max. User Weight: **200kg**

Please contact **REDGUM Brand** if you require further information on this product

## FUNCTIONS



Upholstery



ASK FOR THE BRAND

P: **08 9248 4180**  
E: **orders@for-de.com.au**

[redgumbrand.com.au](http://redgumbrand.com.au)