
COMFORT FROM THE START
JOBST® Elvarex® Soft
DEFY GRAVITY
JOBST
AVAILABLE IN
READY-TO-WEAR
AND
MADE-TO-MEASURE



READY TO WEAR, READY TO GO

JOBST Elvarex Soft is available in a range of **ready-to-wear** styles and sizes, designed to meet a variety of patient needs.

EXPLORE OUR **READY-TO-WEAR** RANGE!



SIMPLIFIED PRODUCT SELECTION WITH OPTIMISED SIZE CHART

The size charts have been optimised to cover more body shapes with fewer sizes. Helping you to save time in clinic and to simplify stock inventory management.
The new sizing system is aligned with the material properties to ensure reliable and consistent medical compression.

SOFTNESS DESIGNED TO SUPPORT CONSISTENT THERAPY

JOBST® Elvarex® Soft is designed to provide effective symptom relief for patients with mild to moderate lymphoedema and lipoedema, where comfort and skin sensitivity are key considerations. The material is designed to provide a soft feeling on the skin and support easy donning, without compromising clinical performance. When compression feels more comfortable, staying consistent with therapy can become easier.

THE TECHNOLOGY BEHIND OUR COMFORT



SOFT TOUCH

The gentle to the touch fabric supports a pleasant feeling on skin using soft yarns. This can enhance wearing comfort, particularly for patients with dry, sensitive or fragile skin.



EASY TO PUT ON

Comfortable daily application and easy to put on due to the combination of knitting technology and soft yarns.



SOFTFIT TECHNOLOGY

The JOBST SoftFit top band helps reduce digging and sliding, offering a more secure and comfortable fit throughout the day. Its finely knitted-in silicone yarn delivers gentle grip with less pressure compared to traditional top bands, helping the garment stay in place without compromising comfort.[2]

1.This is Bhavika's story and not typical of everyone with lymphoedema.
2.Optional for knee high garments and arm sleeves

WHEN TO CHOOSE JOBST ELVAREX SOFT

LYMPHOEDEMA MILD	READY-TO-WEAR
LYMPHOEDEMA MODERATE	MADE-TO-MEASURE
LIPOEDEMA	READY-TO-WEAR OR MADE-TO-MEASURE DEPENDING ON PATIENT NEEDS
ELDERLY OR PEOPLE WITH SENSITIVE SKIN	READY-TO-WEAR OR MADE-TO-MEASURE DEPENDING ON PATIENT NEEDS

APPROVED BY BOTH PATIENTS AND PROFESSIONALS

82% of patients report **high satisfaction** for JOBST Elvarex Soft[1]

TOP RATED BY PROFESSIONALS[1]:

- Soft / lightweight material
- Suitable for sensitive skin
- Comfortable for daily activities
- Easy to apply

1. 3FMR Research: JOBST Lymphology LIVE Segment Portfolio Re-Alignment Qualitative study. 2023

EXPLORE THE READY-TO-WEAR ASSORTMENT

COMPRESSION LEVELS

- 15-21 mmHg / CCL1
- 23-32 mmHg / CCL2
- 34-46 mmHg / CCL3*
  *CCL3 knee high & thigh high only

COLOURS



BLACK

BEIGE

TOP BAND OPTIONS*



SOFTFIT

DOTTED BAND

*Also available with a silicone-free finish.


GARMENT CARE
JOBST ELVAREX SOFT



A laundry bag provides additional protection for the knitted fabric.

*Tumble dryer compatibility proofed

STYLES



KNEE HIGH (AD)

THIGH HIGH (AG)

BERMUDA (ET)

ARMSLEEVE (CG)

GLOVE WITH THUMB AND FINGERS (AC1)

GAUNTLET WITH THUMB (AC1)

Length options
- Petite
- Regular

Toe options
- Open Toe
- Closed Toe

Top band options
- Knitted welt (no silicone)
- 5 m SoftFit silicone band

Length options
- Petite
- Regular

Toe options
- Open Toe
- Closed Toe

Top band options
- 5cm silicone dotted band

Waistband options
- Knitted waistband (no silicone)

Special options
- Female
- Male


Length options
- Petite
- Regular

Top band options
- Knitted welt (no silicone)
- 2.5cm SoftFit silicone band

READY-TO-WEAR WITH **OPTIMISED SIZING**

MEASURING TIPS
FOR ACCURATE SIZING

In order to measure for an optimal fit circumferential measurements should be carried out under different levels of tension. The tension can range from lay-on measurement (0), to pull until 1st resistance (+), to pull until 2nd resistance (++).




All measuring points are relevant for MtM and **the circled measuring points are relevant for RtW.**

FINDING THE RIGHT SIZE
FOR READY-TO-WEAR

With JOBST Elvarex Soft, the sizing has been optimised to cover more body shapes with fewer sizes. This can make finding the right size for your patients more convenient by simplifying product selection.






KNEE HIGH (AD)

Size	1	2	2 plus	3	3 plus	4	4 plus	5
cD (cm)	27 - 34	29 - 37	32 - 40	34 - 43	36 - 46	38 - 49	40 - 51	43 - 56
cC (cm)	28 - 35	30 - 38	33 - 41	35 - 44	37 - 47	39 - 50	41 - 54	44 - 58
cB (cm)	18 - 20	20 - 23	20 - 23	23 - 26	23 - 26	26 - 30	26 - 30	30 - 34


LD Petite: 35 - 40 cm LD Regular: 41 - 46 cm



THIGH HIGH (AG)


Size	1	2	2 plus	3	3 plus	4	4 plus	5
cG (cm)	44 - 57	46 - 60	49 - 63	50 - 66	55 - 71	52 - 70	57 - 75	65 - 82
cD (cm)	27 - 34	29 - 37	32 - 40	34 - 43	36 - 46	38 - 49	40 - 51	43 - 56
cC (cm)	28 - 35	30 - 38	33 - 41	35 - 44	37 - 47	39 - 50	41 - 54	44 - 58
cB (cm)	18 - 20	20 - 23	20 - 23	23 - 26	23 - 26	26 - 30	26 - 30	30 - 34

LG Petite: 68 - 76 cm LG Regular: 77 - 85 cm




BERMUDA (ET)

Size	1	2	3	4	5
cT (cm)	70 - 85	80 - 95	85 - 105	90 - 110	100 - 120
cH (cm)	85 - 100	90 - 105	95 - 115	105 - 125	115 - 135
cG (cm)	44 - 57	46 - 60	50 - 66	52 - 70	65 - 82



GLOVE & GAUNTLET (AC1)

Size	Palm at base of fingers (cA)	Wrist (cC)
1	16.6 - 18 cm	15.6 - 16.7 cm
2	16.6 - 18 cm	16.7 - 18.1 cm
3	18 - 19.5 cm	15.9 - 17 cm
4	18 - 19.5 cm	17 - 18.5 cm
5	18 - 19.5 cm	18.5 - 20.4 cm
6	19.5 - 21.5 cm	17.2 - 18.7 cm
7	19.5 - 21.5 cm	18.7 - 20.6 cm
8	21.5 - 23.1 cm	18.3 - 20.1 cm



ARMSLEEVE (CG)

Size	1	2	3	4	5
cG (cm)	26 - 30	28 - 34	32 - 38	36 - 44	40 - 48
cE (cm)	22 - 26	24 - 28	26 - 32	29 - 35	32 - 40
cC (cm)	14 - 16	16 - 18	18 - 20	20 - 22	22 - 24

LCG Petite: <44 cm LCG Regular: >44 cm

SUPPORT THAT FITS JUST RIGHT

JOBST Elvarex Soft **made-to-measure** is developed to offer a precise fit tailored to each patient's body shape and clinical needs. It is designed to combine reliable compression with easy donning and softness and to perform consistently – even when needs are more complex. This helps support the long-term management of lymphoedema and lipoedema, with confidence for both patients and professionals.



MADE-TO-MEASURE

ASSORTMENT DETAILS MADE-TO-MEASURE

COMPRESSION LEVELS

- 18-21 mmHg / CCL1
- 23-32 mmHg / CCL2
- 34 - 46 mmHg / CCL3
  Lower extremities only

COLOURS

BLACK, BEIGE, MACCHIATO, CRANBERRY, GREY, DARK BLUE, PINE GREEN, DARK BROWN, RUBY RED, YELLOW

TOP BAND OPTIONS*


SOFTFIT

DOTTED BAND

MICRO DOTTED BAND

*Also available with a silicone-free finish.

STYLES



KNEE HIGH (AD) | THIGH HIGH (AG) | PANTY (AT) | CHAP (AG-T) | PANTY SLIP FORM (AT) | CAPRI TIGHTS (BT) | CAPRI TIGHTS (B1T) | CAPRI (CT) | ARMSLEEVE (CG1)

SOFTNESS THAT SUPPORT RESULTS

Delivers **safe and reliable gradient compression** and is designed for a pleasant feeling on skin

Easy to put on to encourage adherence

Supports **simplified product selection** with optimised sizing for ready-to-wear

EXPLORE THE READY-TO-WEAR ASSORTMENT

COMPRESSION LEVELS

18-21 mmHg / CCL1

23-32 mmHg / CCL2

34 - 46 mmHg / CCL3
Lower extremities only

KNEE HIGH (AD)

Length	Top Band	Toe Style	Colour	Size 1	Size 2	Size 2+	Size 3	Size 3+	Size 4	Size 4+	Size 5
Regular	Standard	Closed Toe	Beige	79546-00	79546-01	79546-05	79546-02	79546-06	79546-03	79546-07	79546-04
Regular	Standard	Closed Toe	Black	79546-08	79546-09	79546-13	79546-10	79546-14	79546-11	79546-15	79546-12
Regular	Standard	Open Toe	Beige	79546-32	79546-33	79546-37	79546-34	79546-38	79546-35	79546-39	79546-36
Regular	Standard	Open Toe	Black	79546-40	79546-41	79546-45	79546-42	79546-46	79546-43	79546-47	79546-44
Regular	SoftFit	Closed Toe	Beige	79546-16	79546-17	79546-21	79546-18	79546-22	79546-19	79546-23	79546-20
Regular	SoftFit	Closed Toe	Black	79546-24	79546-25	79546-29	79546-26	79546-30	79546-27	79546-31	79546-28
Regular	SoftFit	Open Toe	Beige	79546-48	79546-49	79546-53	79546-50	79546-54	79546-51	79546-55	79546-52
Regular	SoftFit	Open Toe	Black	79546-56	79546-57	79546-61	79546-58	79546-62	79546-59	79546-63	79546-60
Petite	Standard	Closed Toe	Beige	73605-00	73605-01	73605-05	73605-02	73605-06	73605-03	73605-07	73605-04
Petite	Standard	Closed Toe	Black	73605-08	73605-09	73605-13	73605-10	73605-14	73605-11	73605-15	73605-12
Petite	Standard	Open Toe	Beige	73605-32	73605-33	73605-37	73605-34	73605-38	73605-35	73605-39	73605-36
Petite	Standard	Open Toe	Black	73605-40	73605-41	73605-45	73605-42	73605-46	73605-43	73605-47	73605-44
Petite	SoftFit	Closed Toe	Beige	73605-16	73605-17	73605-21	73605-18	73605-22	73605-19	73605-23	73605-20
Petite	SoftFit	Closed Toe	Black	73605-24	73605-25	73605-29	73605-26	73605-30	73605-27	73605-31	73605-28
Petite	SoftFit	Open Toe	Beige	73605-48	73605-49	73605-53	73605-50	73605-54	73605-51	73605-55	73605-52
Petite	SoftFit	Open Toe	Black	73605-56	73605-57	73605-61	73605-58	73605-62	73605-59	73605-63	73605-60
Regular	Standard	Closed Toe	Beige	79559-24	79559-25	79559-29	79559-26	79559-30	79559-27	79559-31	79559-28
Regular	Standard	Closed Toe	Black	79559-32	79559-33	79559-37	79559-34	79559-38	79559-35	79559-39	79559-36
Regular	Standard	Open Toe	Beige	79559-56	79559-57	79559-61	79559-58	79559-62	79559-59	79559-63	79559-60
Regular	Standard	Open Toe	Black	79559-64	79559-65	79559-69	79559-66	79559-70	79559-67	79559-71	79559-68
Regular	SoftFit	Closed Toe	Beige	79559-40	79559-41	79559-45	79559-42	79559-46	79559-43	79559-47	79559-44
Regular	SoftFit	Closed Toe	Black	79559-48	79559-49	79559-53	79559-50	79559-54	79559-51	79559-55	79559-52
Regular	SoftFit	Open Toe	Beige	79559-72	79559-73	79559-77	79559-74	79559-78	79559-75	79559-79	79559-76
Regular	SoftFit	Open Toe	Black	79559-80	79559-81	79559-85	79559-82	79559-86	79559-83	79559-87	79559-84
Petite	Standard	Closed Toe	Beige	79558-24	79558-25	79558-29	79558-26	79558-30	79558-27	79558-31	79558-28
Petite	Standard	Closed Toe	Black	79558-32	79558-33	79558-37	79558-34	79558-38	79558-35	79558-39	79558-36
Petite	Standard	Open Toe	Beige	79558-56	79558-57	79558-61	79558-58	79558-62	79558-59	79558-63	79558-60
Petite	Standard	Open Toe	Black	79558-64	79558-65	79558-69	79558-66	79558-70	79558-67	79558-71	79558-68
Petite	SoftFit	Closed Toe	Beige	79558-40	79558-41	79558-45	79558-42	79558-46	79558-43	79558-47	79558-44
Petite	SoftFit	Closed Toe	Black	79558-48	79558-49	79558-53	79558-50	79558-54	79558-51	79558-55	79558-52
Petite	SoftFit	Open Toe	Beige	79558-72	79558-73	79558-77	79558-74	79558-78	79558-75	79558-79	79558-76
Petite	SoftFit	Open Toe	Black	79558-80	79558-81	79558-85	79558-82	79558-86	79558-83	79558-87	79558-84
Regular	Standard	Closed Toe	Beige	73552-24	73552-25	73552-29	73552-26	73552-30	73552-27	73552-31	73552-28
Regular	Standard	Closed Toe	Black	73552-32	73552-33	73552-37	73552-34	73552-38	73552-35	73552-39	73552-36
Regular	Standard	Open Toe	Beige	73552-56	73552-57	73552-61	73552-58	73552-62	73552-59	73552-63	73552-60
Regular	Standard	Open Toe	Black	73552-64	73552-65	73552-69	73552-66	73552-70	73552-67	73552-71	73552-68
Regular	SoftFit	Closed Toe	Beige	73552-40	73552-41	73552-45	73552-42	73552-46	73552-43	73552-47	73552-44
Regular	SoftFit	Closed Toe	Black	73552-48	73552-49	73552-53	73552-50	73552-54	73552-51	73552-55	73552-52
Regular	SoftFit	Open Toe	Beige	73552-72	73552-73	73552-77	73552-74	73552-78	73552-75	73552-79	73552-76
Regular	SoftFit	Open Toe	Black	73552-80	73552-81	73552-85	73552-82	73552-86	73552-83	73552-87	73552-84
Petite	Standard	Closed Toe	Beige	73579-24	73579-25	73579-29	73579-26	73579-30	73579-27	73579-31	73579-28
Petite	Standard	Closed Toe	Black	73579-32	73579-33	73579-37	73579-34	73579-38	73579-35	73579-39	73579-36
Petite	Standard	Open Toe	Beige	73579-56	73579-57	73579-61	73579-58	73579-62	73579-59	73579-63	73579-60
Petite	Standard	Open Toe	Black	73579-64	73579-65	73579-69	73579-66	73579-70	73579-67	73579-71	73579-68
Petite	SoftFit	Closed Toe	Beige	73579-40	73579-41	73579-45	73579-42	73579-46	73579-43	73579-47	73579-44
Petite	SoftFit	Closed Toe	Black	73579-48	73579-49	73579-53	73579-50	73579-54	73579-51	73579-55	73579-52
Petite	SoftFit	Open Toe	Beige	73579-72	73579-73	73579-77	73579-74	73579-78	73579-75	73579-79	73579-76
Petite	SoftFit	Open Toe	Black	73579-80	73579-81	73579-85	73579-82	73579-86	73579-83	73579-87	73579-84

JOBST Elvarex Soft RTW knee high garments are supplied as a pair.

THIGH HIGH (AG)

Length	Toe Style	Colour	Size 1	Size 2	Size 2+	Size 3	Size 3+	Size 4	Size 4+	Size 5
Regular	Closed Toe	Beige	73160-00	73160-01	73160-05	73160-02	73160-06	73160-03	73160-07	73160-04
Regular	Closed Toe	Black	73160-08	73160-09	73160-13	73160-10	73160-14	73160-11	73160-15	73160-12
Regular	Open Toe	Beige	73160-16	73160-17	73160-21	73160-18	73160-22	73160-19	73160-23	73160-20
Regular	Open Toe	Black	73160-24	73160-25	73160-29	73160-26	73160-30	73160-27	73160-31	73160-28
Petite	Closed Toe	Beige	73164-00	73164-01	73164-05	73164-02	73164-06	73164-03	73164-07	73164-04
Petite	Closed Toe	Black	73164-08	73164-09	73164-13	73164-10	73164-14	73164-11	73164-15	73164-12
Petite	Open Toe	Beige	73164-16	73164-17	73164-21	73164-18	73164-22	73164-19	73164-23	73164-20
Petite	Open Toe	Black	73164-24	73164-25	73164-29	73164-26	73164-30	73164-27	73164-31	73164-28
Regular	Closed Toe	Beige	73161-00	73161-01	73161-05	73161-02	73161-06	73161-03	73161-07	73161-04
Regular	Closed Toe	Black	73161-08	73161-09	73161-13	73161-10	73161-14	73161-11	73161-15	73161-12
Regular	Open Toe	Beige	73161-16	73161-17	73161-21	73161-18	73161-22	73161-19	73161-23	73161-20
Regular	Open Toe	Black	73161-24	73161-25	73161-29	73161-26	73161-30	73161-27	73161-31	73161-28
Petite	Closed Toe	Beige	73165-00	73165-01	73165-05	73165-02	73165-06	73165-03	73165-07	73165-04
Petite	Closed Toe	Black	73165-08	73165-09	73165-13	73165-10	73165-14	73165-11	73165-15	73165-12
Petite	Open Toe	Beige	73165-16	73165-17	73165-21	73165-18	73165-22	73165-19	73165-23	73165-20
Petite	Open Toe	Black	73165-24	73165-25	73165-29	73165-26	73165-30	73165-27	73165-31	73165-28
Regular	Closed Toe	Beige	73162-00	73162-01	73162-05	73162-02	73162-06	73162-03	73162-07	73162-04
Regular	Closed Toe	Black	73162-08	73162-09	73162-13	73162-10	73162-14	73162-11	73162-15	73162-12
Regular	Open Toe	Beige	73162-16	73162-17	73162-21	73162-18	73162-22	73162-19	73162-23	73162-20
Regular	Open Toe	Black	73162-24	73162-25	73162-29	73162-26	73162-30	73162-27	73162-31	73162-28
Petite	Closed Toe	Beige	73298-00	73298-01	73298-05	73298-02	73298-06	73298-03	73298-07	73298-04
Petite	Closed Toe	Black	73298-08	73298-09	73298-13	73298-10	73298-14	73298-11	73298-15	73298-12
Petite	Open Toe	Beige	73298-16	73298-17	73298-21	73298-18	73298-22	73298-19	73298-23	73298-20
Petite	Open Toe	Black	73298-24	73298-25	73298-29	73298-26	73298-30	73298-27	73298-31	73298-28

JOBST Elvarex Soft RTW thigh high garments are supplied as a pair.

BERMUDA (ET)

Gender	Colour	Size 1	Size 2	Size 3	Size 4	Size 5
Female	Beige	73318-00	73318-01	73318-02	73318-03	73318-04
Female	Black	73318-05	73318-06	73318-07	73318-08	73318-09
Male	Beige	73318-10	73318-11	73318-12	73318-13	73318-14
Male	Black	73318-15	73318-16	73318-17	73318-18	73318-19
Female	Beige	73616-00	73616-01	73616-02	73616-03	73616-04
Female	Black	73616-05	73616-06	73616-07	73616-08	73616-09
Male	Beige	73616-10	73616-11	73616-12	73616-13	73616-14
Male	Black	73616-15	73616-16	73616-17	73616-18	73616-19

JOBST Elvarex Soft RTW bermuda garments are supplied as a single unit.











ARMSLEEVE (CG)

Length	Band	Colour	Size 1	Size 2	Size 3	Size 4	Size 5
Regular	Standard	Beige	73619-00	73619-01	73619-02	73619-03	73619-04
Regular	Standard	Black	73619-05	73619-06	73619-07	73619-08	73619-09
Regular	SoftFit	Beige	73619-10	73619-11	73619-12	73619-13	73619-14
Regular	SoftFit	Black	73619-15	73619-16	73619-17	73619-18	73619-19
Petite	Standard	Beige	78101-02	78101-03	78101-04	78101-05	78101-06
Petite	Standard	Black	78101-07	78101-08	78101-09	78101-10	78101-11
Petite	SoftFit	Beige	78101-12	78101-13	78101-14	78101-15	78101-16
Petite	SoftFit	Black	78101-17	78101-18	78101-19	78101-20	78101-21
Regular	Standard	Beige	78100-02	78100-03	78100-04	78100-05	78100-06
Regular	Standard	Black	78100-07	78100-08	78100-09	78100-10	78100-11
Regular	SoftFit	Beige	78100-12	78100-13	78100-14	78100-15	78100-16
Regular	SoftFit	Black	78100-17	78100-18	78100-19	78100-20	78100-21
Petite	Standard	Beige	78102-03	78102-04	78102-05	78102-06	78102-07
Petite	Standard	Black	78102-08	78102-09	78102-10	78102-11	78102-12
Petite	SoftFit	Beige	78102-13	78102-14	78102-15	78102-16	78102-17
Petite	SoftFit	Black	78102-18	78102-19	78102-20	78102-21	78102-22

JOBST Elvarex Soft RTW armsleeve garments are supplied as a single unit.

GLOVE & GAUNTLET (AC1)

Style	Colour	Size 1	Size 2	Size 3	Size 4	Size 5	Size 6	Size 7	Size 8
Glove	Beige	73156-00	73156-01	73156-02	73156-03	73156-04	73156-05	73156-06	73156-07
Glove	Black	73156-08	73156-09	73156-10	73156-11	73156-12	73156-13	73156-14	73156-15
Gauntlet	Beige	73156-16	73156-17	73156-18	73156-19	73156-20	73156-21	73156-22	73156-23
Gauntlet	Black	73156-24	73156-25	73156-26	73156-27	73156-28	73156-29	73156-30	73156-31
Glove	Beige	73157-00	73157-01	73157-02	73157-03	73157-04	73157-05	73157-06	73157-07
Glove	Black	73157-08	73157-09	73157-10	73157-11	73157-12	73157-13	73157-14	73157-15
Gauntlet	Beige	73157-16	73157-17	73157-18	73157-19	73157-20	73157-21	73157-22	73157-23
Gauntlet	Black	73157-24	73157-25	73157-26	73157-27	73157-28	73157-29	73157-30	73157-31

JOBST Elvarex Soft RTW glove & gauntlet garments are supplied as a single unit.








JOBST® Lymphoedema compression garments

FIND THE RIGHT PRODUCT

ASSESSMENT OF CONDITION

How would you describe your patient's oedema?

- **Mild**
  - What is the level of limb distortion and skinfolds?
    - None to minor
      - What is your patient's tissue like?
        - Soft
        - Firm
    - Minor to moderate → What is your patient's tissue like?
- **Moderate** → What is the level of limb distortion and skinfolds?
- **Severe** → What is the level of limb distortion and skinfolds?
  - Pronounced → What is your patient's tissue like?
    - Soft
    - Very soft to hard/fibrotic/firm

CHALLENGES AND NEEDS

- Soft → Do funding requirements include a ready-to-wear solution, or is accessibility key?[1]
  - Yes → Patient is struggling with donning and/or discomfort due to sensitive skin.
  - No → Patient is struggling with donning and/or discomfort due to sensitive skin.
- Firm → Patient is experiencing discomfort due to feeling restricted while wearing a compression garment. → Patient experiences discomfort due to feeling hot and sweaty.
- Soft → Patient is struggling with donning and/or discomfort due to sensitive skin.
- Very soft to hard/fibrotic/firm → Patient is having difficulty managing their oedema throughout the day and experiences discomfort due to feeling hot. → Difficulty finding the right fit for patients with more complex body shapes.

JOBST® Elvarex® Soft ready-to-wear



Ready-made sizes[1]
Suitable for a variety of needs and with high accessibility.

Easy to put on
Due to the combination of knitting technology and soft yarns.

Soft touch
The gentle to the touch fabric supports a pleasant feeling on skin using soft yarns.

Wall stability of fabric (firm support)

JOBST® Confidence made-to-measure



Freedom of movement
3D ContourFit technology allows the garment to follow the true body shape.

Keeps the skin dry
Advanced moisture management technology.

Wall stability of fabric (firm support)

JOBST® Elvarex® Soft made-to-measure



Easy to put on
Due to the combination of knitting technology and soft yarns.

Soft touch
The gentle to the touch fabric supports a pleasant feeling on skin using soft yarns.

Wall stability of fabric (firm support)

JOBST® Elvarex® made-to-measure



Reliable compression
High containment throughout the day, while helping to bridge skinfolds.


High breathability
Made with an open loop structure.

Wall stability of fabric (firm support)

Developed by Essity based on technical product features.
Individual patient needs and condition must be taken into consideration.

1. RtW, if given sizes do not fit MtM is needed
2. Based on International Society of Lymphology Executive Committee. The Diagnosis and Treatment of Peripheral Lymphedema: 2023 Consensus Document of the International Society of Lymphology. Lymphology. 2023;56(4):133-151

Examples of lymphoedema staging[2]



Stage 0 No swelling | Stage 1 Mild swelling | Stage 2 Moderate swelling | Stage 3 Severe swelling

essity

JOBST®, an Essity brand

ORDER JOBST ELVAREX SOFT

CONTACT ESSITY CLINICAL CONCIERGE

VISIT JOBST ESHOP CUSTOM GARMENT QUOTES & ORDERS

DOWNLOAD CUSTOM GARMENT ORDER FORMS



essity

JOBST®, an Essity brand

Always read the label and follow the directions for use.

ESS1402
07/2026

Essity Health & Medical
PO Box 337, Mount Waverley, Victoria, Australia 3149
medical.essity.com.au
jobst.com.au

Customer Service AU
T 1300 276 633 F 1300 998 830
Customer Service NZ
T 0508 276 111 F 0508 998 830